### EUROPEAN SYNCHROTRON RADIATION FACILITY

INSTALLATION EUROPEENNE DE RAYONNEMENT SYNCHROTRON



## **Experiment Report Form**

The double page inside this form is to be filled in by all users or groups of users who have had access to beam time for measurements at the ESRF.

Once completed, the report should be submitted electronically to the User Office via the User Portal: <a href="https://wwws.esrf.fr/misapps/SMISWebClient/protected/welcome.do">https://wwws.esrf.fr/misapps/SMISWebClient/protected/welcome.do</a>

#### **Deadlines for submission of Experimental Reports**

Experimental reports must be submitted within the period of 3 months after the end of the experiment.

#### **Experiment Report supporting a new proposal ("relevant report")**

If you are submitting a proposal for a new project, or to continue a project for which you have previously been allocated beam time, you must submit a report on each of your previous measurement(s):

- even on those carried out close to the proposal submission deadline (it can be a "preliminary report"),
- even for experiments whose scientific area is different form the scientific area of the new proposal,
- carried out on CRG beamlines.

You must then register the report(s) as "relevant report(s)" in the new application form for beam time.

#### Deadlines for submitting a report supporting a new proposal

- > 1<sup>st</sup> March Proposal Round 5<sup>th</sup> March
- > 10<sup>th</sup> September Proposal Round 13<sup>th</sup> September

The Review Committees reserve the right to reject new proposals from groups who have not reported on the use of beam time allocated previously.

#### Reports on experiments relating to long term projects

Proposers awarded beam time for a long term project are required to submit an interim report at the end of each year, irrespective of the number of shifts of beam time they have used.

#### **Published papers**

All users must give proper credit to ESRF staff members and proper mention to ESRF facilities which were essential for the results described in any ensuing publication. Further, they are obliged to send to the Joint ESRF/ ILL library the complete reference and the abstract of all papers appearing in print, and resulting from the use of the ESRF.

Should you wish to make more general comments on the experiment, please note them on the User Evaluation Form, and send both the Report and the Evaluation Form to the User Office.

#### **Instructions for preparing your Report**

- fill in a separate form for <u>each project</u> or series of measurements.
- type your report in English.
- include the experiment number to which the report refers.
- make sure that the text, tables and figures fit into the space available.
- if your work is published or is in press, you may prefer to paste in the abstract, and add full reference details. If the abstract is in a language other than English, please include an English translation.

<b>ES</b>	R <del>F</del>

Experiment title:	Experiment
X-ray nano-imaging of mitochondria inside neurona	number:
	MD900

LOM		
Beamline:	Date of experiment:	Date of report:
ID16A	from: 12 Apr 2015 to: 17 Apr 2015 from: 7 May 2015 to: 9 May 2015 from: 16 Nov 2015 to: 17 Nov 2015	24/2/21
Shifts:	Local contact(s): Sylvain Bohic	Received at ESRF:

Names and affiliations of applicants (\* indicates experimentalists):

L. Lemelle<sup>1</sup>, A. Simionovici<sup>2,†\*</sup>, P. Colin<sup>3</sup>, G. Knott<sup>4</sup>, S. Bohic<sup>5,6</sup>, P. Cloetens<sup>6</sup>, B.L. Schneider<sup>3\*</sup>

laurence. lemelle@ens-lyon.fr

cloetens@esrf.fr

#### Report: This experimental work was published in

Lemelle, L., Simionovici, A., Colin, P. *et al.* Nano-imaging trace elements at organelle levels in *substantia nigra* overexpressing  $\alpha$ -synuclein to model Parkinson's disease. *Commun Biol* **3**, 364 (2020). https://doi.org/10.1038/s42003-020-1084-0

# Nano-imaging trace elements at organelle levels in *substantia nigra* overexpressing $\alpha$ -synuclein to model Parkinson's disease

Sub-cellular trace element quantifications of nano-heterogeneities in brain tissues offer unprecedented ways to explore at elemental level the interplay between cellular compartments in neurodegenerative pathologies. We designed a quasi-correlative method for analytical nanoimaging of the *substantia nigra*, based on transmission electron microscopy and synchrotron X-ray fluorescence. It combines ultrastructural identifications of cellular compartments and trace element nanoimaging near detection limits, for increased signal-to-noise ratios. Elemental composition of different organelles is compared to cytoplasmic and nuclear compartments in dopaminergic neurons of rat *substantia nigra*. They exhibit 150–460 ppm of Fe, with P/Zn/Fe-rich nucleoli in a P/S-depleted nuclear matrix and Ca-rich rough endoplasmic reticula. Cytoplasm analysis displays sub-micron Fe/S-rich granules, including lipofuscin. Following AAV-mediated overexpression of α-synuclein protein associated with Parkinson's disease, these granules shift towards higher Fe concentrations. This effect advocates for metal (Fe) dyshomeostasis in discrete cytoplasmic regions, illustrating the use of this method to explore neuronal dysfunction in brain diseases.

<sup>&</sup>lt;sup>1</sup>LGL-TPE, ENS de Lyon, Univ. de Lyon, CNRS, 69342 Lyon, France.

<sup>&</sup>lt;sup>2</sup>ISTerre, Univ. Grenoble Alpes, Univ. Savoie Mont Blanc, CNRS, IRD, IFSTTAR, CS 40700, 38058 Grenoble, France

<sup>&</sup>lt;sup>†</sup>Institut Universitaire de France (IUF). <u>alexandre.simionovici@univ-grenoble-alpes.fr</u>

<sup>&</sup>lt;sup>3</sup>Brain Mind Institute, Ecole Polytechnique Fédérale de Lausanne (EPFL), 1015 Lausanne, Switzerland. <a href="mailto:philippe.colin@epfl.ch">philippe.colin@epfl.ch</a>, bernard.schneider@epfl.ch

<sup>&</sup>lt;sup>4</sup>Centre of Interdisciplinary Electron Microscopy, Ecole Polytechnique Fédérale de Lausanne (EPFL), 1005 Lausanne, Switzerland. graham.knott@epfl.ch

<sup>&</sup>lt;sup>5</sup>INSERM UA7, Rayonnement Synchrotron pour la Recherche Biomédicale, STROBE, 38043 Grenoble, France. bohic@esrf.fr <sup>6</sup>ESRF - The European Synchrotron, ID16A beamline, 38043 Grenoble Cedex 9, France.